

DRSHD 1080p™

The Med X Change, Inc. DRSHD 1080p™ offers the latest in FULL 1080p video recording with a full featured GUI (Graphical User Interface), all while delivering brilliant HD video quality.

FULL 1080P HD VIDEO RECORDING

Support for Full 1080p Video Recording with DVI input for 1080p, 1080i, 720p NTSC & PAL

CAPTURE HD VIDEO & STILL FROM HD MICROSCOPE/ENDOSCOPE

H.264 .mp4/.mov/.avi, JPEG, BMP, PNG, DICOM

VOICE TO TEXT ANNOTATIONS

Customizable annotations

DICOM SCU INTEGRATION

For Still Images and H.264 HD Video



PRINTING/PLAYBACK/COPY

- Auto Print to Sony Medical Grade Printers
- Playback Video/Stills & Copy to external USB Media
- Still Image Enhancements & Voice to Text Annotations



HIPPA SECURE FILES

- Optional security feature to prevent unauthorized access to software application
- Patient information is only linked to media within the software application database
- Users have the option to remove patient information from media through copy, burn functions and archived cases

Power Requirements: 5V DCA, 25w
Medical Compliance: CE, TUV, CAN/CSA 22.2 No. 601.1-M90
UL60601-1:2003, IEC 60601-1-2:2007, IEC 3rd Edition

MODALITY WORKLIST INTEGRATION

MULTI LANGUAGE SUPPORT ONE TOUCH RECORDING




MEDXSTREAM iOS & ANDROID APP

Pair your phone/tablet to your DRSHD 1080p™ with our MedXStream™ app where you can instantly access your Surgical Cases to Review, Annotate, Playback, and Download. Visit medxstream.com for more information.



STORAGE

- Large Storage Capacity – 2 TB hard drive, approximately 170 hours of HD video
- Record to internal and external USB hard drive



POST PRODUCTION

- Create Surgical Reports arrows and annotations through the MedXStream™ app
- Native Mac & Windows 7, 8 Editing – Edit with Final Cut Pro, Quicktime, iMovie or Adobe Premier Pro
- Easily add text annotations to videos and still images during recording or editing processes

ISO 13485 Certified

©2015 Med X Change, Inc. All Rights Reserved. All trademarks and registered trademarks are the property of their respective owners.



Med X Change, Inc.
medxchange.com

24hr Tech Support: 1 (888) 982-4469
support@medxchange.com